

**EFFICACY DATA for ES64H**  
**BTC<sup>®</sup> 885 Neutral Disinfectant Cleaner-64**  
**(EPA Reg. No. 1839-169)**

**DISINFECTION DATA:**

**Test Method:** AOAC Use Dilution

**Test Conditions:** 10 minute contact time, 5% organic soil load, stainless steel carrier substrates, 400 ppm hard water, 20°C exposure temperature, 2 oz/gal dilution [with the exception of *Xanthomonas axonopodis* pathovar *citri*, dilution ratio is indicated by organism]

**Results:**

<u>Test Organism</u>	<u>Sample</u>	<u>No. of Carriers</u>	
		<u>Exposed</u>	<u>Positive</u>
<i>Staphylococcus aureus</i> (ATCC 6538)	A	60	1
	B	60	1
	C	60	1
<i>Salmonella (choleraesuis) enteric</i> (ATCC 10708)	A	60	0
	B	60	0
	C	60	1
<i>Pseudomonas aeruginosa</i> PRD-10 (ATCC 15442)	A	60	0
	B	60	0
	C	60	1
Ampicillin resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
Bactrim resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
<i>Bordetella bronchiseptica</i> (ATCC 31437)	A	10	0
	B	10	0
Cefazolin resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
Ceftazidime resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
Ceftriaxone resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
Ciprofloxacin resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
Community Associated Methicillin Resistant <i>Staphylococcus aureus</i> (CA-MRSA) (NRS 123, Genotype USA400)	A	10	0
	B	10	0
	C	10	0
Community Associated Methicillin Resistant <i>Staphylococcus aureus</i> (CA-MRSA) (NRS 384, Genotype USA300)	A	10	0
	B	10	0
	C	10	0
<i>Corynebacterium ammoniagenes</i> (ATCC 6871)	A	10	0
	B	10	0
<i>Enterobacter aerogenes</i> (ATCC 13048)	A	10	0
	B	10	0
<i>Enterobacter cloacae</i> (ATCC 23355)	A	10	0
	B	10	0
<i>Enterobacter cloacae</i> (clinical isolate)	A	10	0
	B	10	0

- continued on next page -

**DISINFECTION DATA** (continued)

<u>Test Organism</u>	<u>Sample</u>	<u>No. of Carriers</u>	
		<u>Exposed</u>	<u>Positive</u>
<i>Enterococcus faecalis</i> (ATCC 19433)	A	10	0
	B	10	0
<i>Enterococcus faecalis</i> (clinical isolate)	A	10	0
	B	10	0
<i>Escherichia coli</i> (ATCC 11229)	A	10	0
	B	10	0
<i>Escherichia coli</i> (clinical isolate)	A	10	0
	B	10	0
<i>Escherichia coli</i> O111:H8 (ATCC BAA-184)	A	10	0
	B	10	0
<i>Fusobacterium necrophorum</i> (ATCC 27852)	A	10	0
	B	10	0
Gentamicin resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
<i>Klebsiella pneumoniae</i> subsp. <i>Pneumonia</i> (ATCC 13883)	A	10	0
	B	10	0
<i>Lactobacillus casei</i> subsp. <i>Rhamnosus</i> (ATCC 7469)	A	10	0
	B	10	0
Levofloxacin resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
<i>Listeria monocytogenes</i> (ATCC 35152)	A	10	0
	B	10	0
Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA) (ATCC 33592)	A	10	0
	B	10	0
<i>Pasteurella multocida</i> (ATCC 7707)	A	10	0
	B	10	0
<i>Proteus mirabilis</i> (ATCC 9921)	A	10	0
	B	10	0
<i>Proteus mirabilis</i> (ATCC 25933)	A	10	0
	B	10	0
<i>Proteus vulgaris</i> (ATCC 13315)	A	10	0
	B	10	0
<i>Salmonella (paratyphi B) enterica</i> (ATCC 8759)	A	10	0
	B	10	0
<i>Salmonella (pullorum) enterica</i> (ATCC 9120)	A	10	0
	B	10	0
<i>Salmonella (typhi) enterica</i> (ATCC 6539)	A	10	0
	B	10	0
<i>Salmonella (typhimurium) enterica</i> (ATCC 14028)	A	10	0
	B	10	0
<i>Salmonella (enteritidis) enterica</i> (ATCC 13076)	A	10	0
	B	10	0
<i>Serratia marcescens</i> (ATCC 8100)	A	10	0
	B	10	0
<i>Shigella dysenteriae</i> (ATCC 12180)	A	10	0
	B	10	0
<i>Shigella flexneri</i> Type 2b (ATCC 12022)	A	10	0
	B	10	0

- continued on next page -

DISINFECTION DATA (continued)

<u>Test Organism</u>	<u>Sample</u>	<u>No. of Carriers</u>	
		<u>Exposed</u>	<u>Positive</u>
<i>Shigella sonnei</i> (ATCC 25931)	A	10	0
	B	10	0
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (ATCC 33592)	A	10	0
	B	10	0
<i>Staphylococcus aureus</i> (clinical isolate)	A	10	0
	B	10	0
<i>Staphylococcus epidermidis</i> (ATCC 29641)	A	10	0
	B	10	0
<i>Staphylococcus epidermidis</i> (clinical isolate)	A	10	0
	B	10	0
<i>Streptococcus pyogenes</i> Group A (ATCC 19615)	A	10	0
	B	10	0
<i>Streptococcus pyogenes</i> (clinical-flesh eating strain, BIRD M3)	A	10	0
	B	10	0
Tobramycin resistant <i>Acinetobacter baumannii</i> (Fairfax Hospital CI 02001)	A	10	0
	B	10	0
Vancomycin Resistant <i>Enterococcus faecalis</i> (VRE) (ATCC 51575)	A	10	0
	B	10	0
Vancomycin Intermediate Resistant <i>Staphylococcus aureus</i> (VISA) (HIP 5863)	A	10	0
	B	10	0
<i>Xanthomonas maltophilia</i> (clinical isolate)	A	10	0
	B	10	0
<i>Xanthomonas axonopodis</i> pathovar <i>citri</i> at a dilution ratio of 1:27 (2000 ppm active quat)	A	10	0
	B	10	0

**Conclusion:** Under the conditions of these investigations, BTC® 885 Neutral Disinfectant Cleaner-64 demonstrated **disinfectant** activity against *Staphylococcus aureus*, *Salmonella (choleraesuis) enterica*, *Pseudomonas aeruginosa*, Ampicillin resistant *Acinetobacter baumannii*, Bactrim resistant *Acinetobacter baumannii*, *Bordetella bronchiseptica*, Cefazolin resistant *Acinetobacter baumannii*, Ceftazidime resistant *Acinetobacter baumannii*, Ceftriaxone resistant *Acinetobacter baumannii*, Ciprofloxacin resistant *Acinetobacter baumannii*, Community Associated Methicillin Resistant *Staphylococcus aureus* (CA-MRSA) (NRS 123, Genotype USA400), Community Associated Methicillin Resistant *Staphylococcus aureus* (CA-MRSA) (NRS 384, Genotype USA300), *Corynebacterium ammoniagenes*, *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* O111:H8, *Fusobacterium necrophorum*, Gentamicin resistant *Acinetobacter baumannii*, *Klebsiella pneumoniae* subsp. *pneumoniae*, *Lactobacillus casei* subsp. *rhamnosus*, Levofloxacin resistant *Acinetobacter baumannii*, *Listeria monocytogenes*, Methicillin Resistant *Staphylococcus aureus* (MRSA), *Pasteurella multocida*, *Proteus mirabilis* (ATCC 9921), *Proteus mirabilis* (ATCC 25933), *Proteus vulgaris*, *Salmonella (paratyphi B) enterica*, *Salmonella (pullorum) enterica*, *Salmonella (typhi) enterica*, *Salmonella (typhimurium) enterica*, *Salmonella (enteritidis) enterica*, *Serratia marcescens*, *Shigella dysenteriae*, *Shigella flexneri* Type 2b, *Shigella sonnei*, *Staphylococcus aureus* subsp. *aureus*, *Staphylococcus epidermidis*, *Streptococcus pyogenes* Group A, *Streptococcus pyogenes* (clinical-flesh eating strain, BIRD M3), Tobramycin resistant *Acinetobacter baumannii*, Vancomycin Resistant *Enterococcus faecalis* (VRE), and Vancomycin Intermediate Resistant *Staphylococcus aureus* (VISA), according to criteria established by the U.S. Environmental Protection Agency for registration and labeling of a disinfectant product as a bactericide.

BTC® 885 Neutral Disinfectant Cleaner-64 also demonstrated **disinfectant** activity against the following antibiotic resistant clinical isolates: *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, and *Xanthomonas maltophilia*.

At a dilution ratio of 1:27 (2000 ppm active quat), BTC® 885 Neutral Disinfectant Cleaner-64 demonstrated **disinfectant** activity against *Xanthomonas axonopodis* pathovar *citri* (Citrus Canker Disease).

**VIRUCIDAL DATA:**

**Test Method:**

- \* U.S. E.P.A. Pesticide Assessment Guidelines, Subdivision G: Product Performance, Section 91-2 (f), and Section 91-30, (d), (e), November 1982.
- † Protocols for Testing the Efficacy of Disinfectants against Hepatitis B Virus (HBV) (EPA, Federal Register, Vol. 65, No. 166, 8/25/2000, p. 51828).
- ‡ Protocol for Testing Disinfectants against Hepatitis C Virus using Bovine Viral Diarrhea Virus as approved by the U.S. EPA on August 15, 2002.

**Test Conditions:** 2 oz/gal dilution, 5% organic soil load, 10 minute contact time,  
400 ppm hard water, 21-24°C exposure temperature, sterile glass petri dishes,

**Results:**

Test Organism	Sample		Titer Reduction	
	A	B		
*Avian Influenza A Virus (H3N2) (Avian Reassortant) (ATCC VR-2072)	A	B	≥4.25 log <sub>10</sub>	≥4.25 log <sub>10</sub>
*Avian Influenza Virus, Type A (Turkey/WIS/66) (H9N2)	A	B	≥4.0 log <sub>10</sub>	≥4.0 log <sub>10</sub>
*Bovine Rhinotracheitis, strain LA (ATCC VR-188)	A	B	≥5.0 log <sub>10</sub>	≥5.0 log <sub>10</sub>
‡Bovine Viral Diarrhea Virus (BVDV)	A	B	5.9 log <sub>10</sub>	5.9 log <sub>10</sub>
*Canine Distemper Virus, strain Lederle (ATCC VR-128)	A	B	≥6.25 log <sub>10</sub>	≥6.25 log <sub>10</sub>
*Feline Picornavirus, strain FRV (ATCC VR-649)	A	B	≥4.25 log <sub>10</sub>	≥4.25 log <sub>10</sub>
†Hepatitis B Virus (HBV) (Duck Hepatitis B Virus-DHBV)	A	B	4.5 log <sub>10</sub>	4.7 log <sub>10</sub>
‡Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus-BVDV)	A	B	5.9 log <sub>10</sub>	5.9 log <sub>10</sub>
*Herpes Simplex Type 1 (ATCC VR-260)	A	B	≥5.0 log <sub>10</sub>	≥5.0 log <sub>10</sub>
*Herpes Simplex Type 2 (ATCC VR-734)	A	B	≥6.0 log <sub>10</sub>	≥6.0 log <sub>10</sub>
*Human Coronavirus (ATCC VR-740, strain 229E)	A	B	≥3.0 log <sub>10</sub>	≥3.0 log <sub>10</sub>
*Human Immunodeficiency Virus, HTLV-III <sub>RF</sub> , strain of HIV-1 (associated with AIDS)	A	B	≥3.5 log <sub>10</sub>	≥3.5 log <sub>10</sub>
*Human Immunodeficiency Virus type 2 (HIV-2), strain CBL-20	A	B	≥3.25 log <sub>10</sub>	≥3.25 log <sub>10</sub>
*Influenza A <sub>2</sub> , strain Hong Kong (ATCC VR-544)	A	B	≥4.25 log <sub>10</sub>	≥4.25 log <sub>10</sub>
*Pandemic 2009 H1N1 Influenza A Virus	(Refer to NOTE below.)			
*Paramyxovirus (Mumps) (ATCC VR-1438)	A	B	≥3.0 log <sub>10</sub>	≥3.0 log <sub>10</sub>
*Porcine Respiratory & Reproductive Syndrome Virus (PRRSV), strain NVSL	A	B	≥5.0 log <sub>10</sub>	≥5.0 log <sub>10</sub>
*Pseudorabies, strain Aujeszky (ATCC VR-135)	A	B	≥5.25 log <sub>10</sub>	≥5.25 log <sub>10</sub>
*Rabies Virus (attenuated CDC ERA strain)	A	B	3.0 log <sub>10</sub>	3.0 log <sub>10</sub>
*Rotavirus, strain SA-11 (ATCC VR-899)	A	B	4.5 log <sub>10</sub>	4.5 log <sub>10</sub>
*SARS Associated Coronavirus (ZeptoMetrix)	A	B	3.03 log <sub>10</sub>	3.03 log <sub>10</sub>
*Vaccinia, strain WR (ATCC VR-119)	A	B	≥5.5 log <sub>10</sub>	≥5.5 log <sub>10</sub>

**Conclusion:** Under the conditions of this investigation, BTC® 885 Neutral Disinfectant Cleaner-64 demonstrated **virucidal** activity against Avian Influenza A Virus (H3N2), Avian Influenza Virus Type A (H9N2), Bovine Rhinotracheitis, Bovine Viral Diarrhea Virus (BVDV), Canine Distemper Virus, Feline Picornavirus, Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), Herpes Simplex Type 1, Herpes Simplex Type 2, Human Coronavirus, Human Immunodeficiency Virus (HIV-1), Human Immunodeficiency Virus type 2 (HIV-2), Influenza A<sub>2</sub>, Pandemic 2009 H1N1 Influenza A Virus, Paramyxovirus (Mumps), Porcine Respiratory & Reproductive Syndrome Virus (PRRSV), Pseudorabies, Rabies Virus, Rotavirus, SARS Associated Coronavirus and Vaccinia according to criteria established by the U. S. Environmental Protection Agency for registration and labeling of a disinfectant product as a virucide.

**NOTE:** Per the EPA guidance document dated October 21, 2009, disinfectant products that bear label claims against human, avian, or swine influenza A virus, and have submitted and received approval of efficacy data to support these label claims, may include a label claim against the Pandemic 2009 H1N1 Influenza A Virus.

**FUNGICIDAL DATA:**

**Test Method:** AOAC Fungicidal Activity of Disinfectants

**Test Conditions:** 2 oz/gal dilution, 5% organic soil load, 10 minute contact time,  
 200 ppm hard water, 20°C exposure temperature

<u>Test Organism</u>	<u>Sample</u>	<u>Exposure Time</u> (min.) vs. Growth		
		<u>5</u>	<u>10</u>	<u>15</u>
<i>Trichophyton mentagrophytes</i> (ATCC 9533)	A	+	0	0
	B	+	0	0
<i>Candida albicans</i> (ATCC 10231)	A	0	0	0
	B	0	0	0

**Conclusions:** Under the conditions of this investigation, BTC® 885 Neutral Disinfectant Cleaner-64 demonstrated **fungicidal** activity against *Trichophyton mentagrophytes* and *Candida albicans* according to criteria established by the U. S. Environmental Protection Agency for registration and labeling of a disinfectant product as a fungicide.

**MILDEW FUNGISTATIC DATA:**

**Test Method:** Hard Surface Mildew Fungistatic Test

**Test Organism:** *Aspergillus niger* (ATCC 6275)

**Test Conditions:** 2 oz/gal dilution, 400 ppm hard water, ceramic tile carriers

**Results:**

<u>Sample</u>	<u>No. of Exposed Tiles</u>	<u>No. of Tiles Showing Growth</u>
A	10	0
B	10	0
Control	10	10

**Conclusion:** Under the conditions of this investigation, BTC® 885 Neutral Disinfectant Cleaner-64 demonstrated **fungistatic** activity against *Aspergillus niger* according to criteria established by the U. S. Environmental Protection Agency for registration and labeling of a disinfectant product as a fungistat.

**BACTERICIDAL STABILITY DATA OF USE-SOLUTION:**

**Test Method:** Use Dilution

**Test Conditions:** 2 oz/gal dilution, 5% organic soil load, 10 minute contact time,  
 400 ppm hard water, 20°C exposure temperature, stainless steel carrier substrates

**Storage Conditions:** sealed containers at room temperature

**Results:**

<u>Test Time</u>	<u>Sample</u>	<u>Test Organism</u>	<u>No. of Carriers</u>	
			<u>Exposed</u>	<u>Positive</u>
Zero Time	A	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
	B	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
Week 1	A	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
	B	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
Week 2	A	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
	B	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
Week 3	A	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
	B	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
Week 4	A	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0
	B	<i>Staphylococcus aureus</i> (ATCC 6538)	10	0
		<i>Salmonella (choleraesuis) enterica</i> (ATCC 10708)	10	0
		<i>Pseudomonas aeruginosa</i> (ATCC 15442)	10	0

**Conclusion:** The results of this investigation show that a 2 oz/gal use dilution of BTC® 885 Neutral Disinfectant Cleaner-64 will demonstrate disinfectant efficacy against *Staphylococcus aureus*, *Salmonella (choleraesuis) enterica*, and *Pseudomonas aeruginosa* for up to 4 weeks in accordance with criteria established by the U.S. Environmental Protection Agency for registration and labeling of a disinfectant product as a bactericide.