

## TASKMISER DISPENSER

Dilution proportioning tip chart for the 1.89L/0.50 U.S. Gal. bottle used in the Taskmiser Dispenser

ES PRODUCT	SINGLE FLOW		DUAL FLOW		
	HF DILUTION	LF DILUTION	HF TIP	LF TIP	OFF HF TIP OF BOTTLE + TIP @ LF EDUCATOR
70C - GP Cleaner	1:30	1:20	White	Orange	Purple; 1:64 LF, Olive Green; 1:128 HF + 1:30 LF
71+C - Super H2O2	1:128	1:20	Aqua	Brown	n/a
74C - Spray Et Wipe	1:16	1:7	No Tip	Red	Red; 1:7 LF Et 1:16 HF
75C - HD Cleaner/Degreaser	1:32	1:20	White	Orange	Green; 1:20 LF Et 1:32 HF
77C NFP - Glass Cleaner	1:32	1:32	Blue	Yellow	n/a
51C - Washroom Cleaner	1:32	1:16	No Tip	Red	No Tip; 1:16 LF Et 1:33 HF
56C Mild Acid Washroom Cleaner	1:16	1:5	No Tip	No Tip	Red; 1:7 LF Et 1:16 HF
58C - Odor Eliminator	1:16	1:5	No Tip	No Tip	Red; 1:7 LF Et 1:16 HF
83C - Haze Away	1:128	1:100	Light Purple	Aqua	n/a
84C - Neutral Floor Cleaner	1:267	1:236	Light Grey	Light Orange	Light Purple; 1:236 LF Et 1:267 HF
87C - Carpet Extraction Cleaner with Envirocide	1:64	1:32	Orange	Aqua	n/a
90C - Clean Cut Stripper	1:20	1:13	n/a	Light Blue	n/a
91C - Deep Scrub	1:38	1:32	Light Blue	Aqua	n/a
ES95 Restoration	1:36	1:38	Blue	Aqua	n/a
64C or 64HC- GP Neutral Disinfectant Cleaner	1:64	1:56	Orange	Purple	New Red; 1:63 LF Et 1:64 HF
256HC - Neutral Disinfectant *Hospital Grade	1:256	1:256	Purple	Light Purple	n/a
512C - Sanitizer/ Disinfectant	1:490	1:150	Red-Purple	Light Orange	Light Orange; 1:503 LF Et 1:490 HF
25C - Disinfectant/HD Cleaner/ Odor Eliminator	1:26	1:22	Beige	Green	Light Blue; 1:17 LF Et 1:26 HF
128C - Disinfectant-Sanitizer-Fungicide-Deodorizer	1:110	1:120	Aqua	Pink	n/a
920C - Bio-Active Bathroom Cleaner	1:330	1:370	New Black	Light Orange	n/a
680C - Bio-Active Super Floor Degreaser	1:445	1:220	Pink	Olive	n/a